



## Microlab Technologies - Consumables Kit

**Model:**

**Serial No:**

**Valves** – Should be changed as required (usually indicated by air intake or poor performance)

Valve Block, Single Dispenser, Model: ML501A, ML510B, ML512BP	Ref: <b>35825</b>	<input type="checkbox"/>
Valve Block, Diluter Valve, Model: ML503A, ML530B, ML532BP	Ref: <b>35844</b>	<input type="checkbox"/>
Valve Block, Dispenser Valve, Model: ML504A, ML540B, ML542BP	Ref: <b>35842</b>	<input type="checkbox"/>

**Syringes** – Should be changed every 6 months

### Left or Right Side Diluter/ Dispenser

25µl TLLX Syringe, Ref: <b>80222</b>	<input type="checkbox"/>
50µl TLLX Syringe, Ref: <b>80922</b>	<input type="checkbox"/>
100µl TLLX Syringe, Ref: <b>81022</b>	<input type="checkbox"/>
250µl TLLX Syringe, Ref: <b>81122</b>	<input type="checkbox"/>
500µl TLLX Syringe, Ref: <b>81222</b>	<input type="checkbox"/>
1000µl TLLX Syringe, Ref: <b>81323</b>	<input type="checkbox"/>
1000µl TLL Syringe, Ref: <b>81320</b>	<input type="checkbox"/>
2500µl TLL Syringe, Ref: <b>81420</b>	<input type="checkbox"/>
5000µl TLL Syringe, Ref: <b>81520</b>	<input type="checkbox"/>
10000µl TLL Syringe, Ref: <b>81620</b>	<input type="checkbox"/>
25000µl TLL Syringe, Ref: <b>82521</b>	<input type="checkbox"/>

### Right Side Diluter

25µl DX Syringe, Ref: <b>80226</b>	<input type="checkbox"/>
50µl DX Syringe, Ref: <b>80926</b>	<input type="checkbox"/>
100µl DX Syringe, Ref: <b>81026</b>	<input type="checkbox"/>
250µl DX Syringe, Ref: <b>81126</b>	<input type="checkbox"/>
500µl DX Syringe, Ref: <b>81226</b>	<input type="checkbox"/>
1000µl DX Syringe, Ref: <b>81326</b>	<input type="checkbox"/>

**Salt line syringes** ( for saline solutions that crystallize & scratch traditional PTFE plunger tips.)

250µl TLL-SAL Syringe, Ref: <b>203220</b>	<input type="checkbox"/>	250µl SAL-D, Ref: <b>203450</b>	<input type="checkbox"/>
500µl TLL-SAL Syringe, Ref: <b>203230</b>	<input type="checkbox"/>	500µl SAL-D, Ref: <b>203460</b>	<input type="checkbox"/>
1000µl TLL-SAL Syringe, Ref: <b>203240</b>	<input type="checkbox"/>	1000µl SAL-D, Ref: <b>203470</b>	<input type="checkbox"/>
2500µl TLL-SAL Syringe, Ref: <b>203250</b>	<input type="checkbox"/>		
5000µl TLL-SAL Syringe, Ref: <b>203260</b>	<input type="checkbox"/>		
10000µl TLL-SAL Syringe, Ref: <b>203270</b>	<input type="checkbox"/>		

**Tubing:** Should be changed every 3 months

**Micro** (18 gauge) for syringe sizes below 2500µl

**Fill Left Side,** Ref: **240010**

**Dispense Right Side,** Ref: **240130**

**Macro** (12 gauge) for syringe sizes above 2500µl

**Fill Left Side,** Ref: **240000**

**Dispense Right Side,** Ref: **240360**

### Accessories:

Tubing Clip, Blue	Ref: <b>230010</b>	<input type="checkbox"/>
Sml Parts Kit	Ref: <b>35888</b>	<input type="checkbox"/>
Tubing Grip	Ref: <b>342-9892</b>	<input type="checkbox"/>

Hand Grip Holder	Ref: <b>35783</b>	<input type="checkbox"/>
Hand Actuator	Ref: <b>35529</b>	<input type="checkbox"/>
Dispenser Actuator	Ref: <b>35767</b>	<input type="checkbox"/>



Microlab Technologies trained & approved service agents for Hamilton Microlab® products.

Microlab Technologies Ltd • Orchard Lodge • Royer Close • Hawkwell • Essex SS5 4LR

Telephone: +44(0)1702 208044 • Fax: +44(0)1702 202364 • Email: sales@microlabtech.co.uk • Web: www.microlabtech.co.uk